


ELAST-O-TIES

These latex-free elastic ligature ties are designed to maximize office efficiency while promoting patient compliance. Our features include:

- ◆ A very high resistance to breakage or tearing during placement or treatment, saving valuable chair time and promoting patient confidence.
- ◆ High flexibility and resiliency without stiffness, making the ties easier to place and remove while providing a longer lasting, more effective tie.
- ◆ Our ties are colorfast to improve patient satisfaction and compliance.
- ◆ 84 ties per stick, 12 sticks per pack, for a total of 1,008 ties per pack.



COLORS

GRAY	300-100	MAGENTA	300-117
CLEAR	300-101	KELLY GREEN	300-118
PINK	300-102	MARIGOLD	300-119
NEON GREEN	300-103	SPARKLE GREEN	300-120
PURPLE	300-104	SMOKE GRAY	300-121
RED	300-105	TEAL	300-122
ROYAL BLUE	300-106	FUCHSIA	300-123
BLACK	300-107	MIDNIGHT BLUE	300-124
ORANGE	300-108	GLOW IN DARK	300-125
YELLOW	300-109	PERIWINKLE	300-126
WHITE	300-110	CAROLINA BLUE	300-127
IVORY/TOOTH	300-111	BABY BLUE	300-128
EMERALD	300-112	SPARKLE RED	300-129
TURQUOISE	300-113	SPARKLE BLUE	300-130
LIGHT BLUE	300-114	BABY PINK	300-131
SILVER	300-115	LILAC	300-132
GOLD	300-116	PEARL	300-133
		CORAL	300-134

ROTATION WEDGES

Our elastomeric rotation wedges fit snugly and have knobs to ensure they stay underneath the archwire. The handle of the dispenser is wedge-shaped for quick identification on your chairside table. They are sold in gray and clear. They contain no latex, BPA or plasticizers.



GRAY	300-500
CLEAR	300-501

BRACKET ORGANIZERS

CATALOG # 100-001



PATIENT-SIZE TIE PACKS


BITE SIZE TIES



42 pieces per package,
for a total of
1,008 ties.

BITE-SIZE TIES	COLORS	TEN-TIES
300-200	GRAY	300-300
300-201	CLEAR	300-301
300-202	PINK	300-302
300-203	GREEN	300-303
300-204	PURPLE	300-304
300-205	RED	300-305
300-206	BLUE	300-306
300-207	BLACK	300-307
300-208	ORANGE	300-308
300-209	YELLOW	300-309
300-210	WHITE	300-310
300-211	IVORY	300-311
300-212	EMERALD	300-312
300-213	TURQUOISE	300-313
300-214	LIGHT BLUE	300-314
300-215	SILVER	300-315
300-216	GOLD	300-316
300-217	MAGENTA	300-317
300-218	KELLY GREEN	300-318
300-219	MARIGOLD	300-319
300-220	SPARKLE GREEN	300-320
300-221	SMOKE GRAY	300-321
300-222	TEAL	300-322
300-223	FUCHSIA	300-323
300-224	MIDNIGHT BLUE	300-324
300-225	GLOW IN DARK	300-325
300-226	PERIWINKLE	300-326
300-227	CAROLINA BLUE	300-327
300-228	BABY BLUE	300-328
300-229	SPARKLE RED	300-329
300-230	SPARKLE BLUE	300-330
300-231	BABY PINK	300-331
300-232	LILAC	300-332
300-233	PEARL	300-333
300-234	CORAL	300-334

TEN - TIES



100 pieces per package, for
a total of 1,000 ties.
Contains no latex!

ORTHODONTIC ELASTICS (LATEX and NON-LATEX)

All catalog numbers are given for latex.

For non-latex, add "NL" to number when ordering.



PATIENT-SIZED QUANTITIES

50 bags of 100 elastics/bag (5,000 elastics)

Diameter	1/8"	3/16"	1/4"	5/16"	3/8"	COLOR
Bag Color	Purple	Green	Red	Blue	Orange	
2.5oz. Light	811-000	811-003	811-006	811-009	811-012	AMBER
4.5oz. Med	811-001	811-004	811-007	811-010	811-013	
6.5oz Heavy	811-002	811-005	811-008	811-011	811-014	
2.5oz. Light	811-175	811-176	811-177	811-178	811-179	NEON
4.5oz. Med	811-015	811-017	811-019	811-021	811-023	
6.5oz Heavy	811-016	811-018	811-020	811-022	811-024	

BULK QUANTITIES

Our standard elastics bagged bulk for your special dispensing requirements.

Bags of 10,000

Diameter	1/8"	3/16"	1/4"	5/16"	3/8"	COLOR
2.5oz. Light	811-025	811-028	811-031	811-034	811-037	
4.5oz. Med	811-026	811-029	811-032	811-035	811-038	
6.5oz Heavy	811-027	811-030	811-033	811-036	811-039	
2.5oz. Light	811-170	811-171	811-172	811-173	811-174	NEON
4.5oz. Med	811-040	811-042	811-044	811-046	811-048	
6.5oz Heavy	811-041	811-043	811-045	811-047	811-049	

Bags of 100,000

Diameter	1/8"	3/16"	1/4"	5/16"	3/8"	COLOR
2.5oz. Light	811-200	811-203	811-206	811-209	811-212	
4.5oz. Med	811-201	811-204	811-207	811-210	811-213	
6.5oz Heavy	811-202	811-205	811-208	811-211	811-214	
2.5oz. Light	811-230	811-231	811-232	811-233	811-234	NEON
4.5oz. Med	811-220	811-222	811-224	811-226	811-228	
6.5oz Heavy	811-221	811-223	811-225	811-227	811-229	

EXTRA-ORAL ELASTICS, Bags of 5,000

Diameter	1/4"	5/16"	3/8"	1/2"	COLOR
16oz.	811-125	811-126	811-127	811-128	
8oz.	811-056	811-057	811-058	811-059	

SEP-O-LOOPS

- * .180" in Diameter
- * Stretchable Material
- * Excellent Memory
- * Ideal for Posterior Separators
- * Will Stretch out to Maximum Limit using Separator Plier
- * Will Return Creating Spacing Between Teeth

SEPO IIs (SEP-O-LOOPS with knobs)




Molded *ROUND* edges—unlike cut tubing separators and made only for orthodontists

Molded 80 modules per unit—No fumbling

Knobs help keep subgingival loss from occurring hence fewer nuisance radiographs

SEPO-II separators come in two thicknesses for different degrees of tooth spacing. Use the Purple separators for teeth requiring little movement before banding. Our Light Blue separators are for instances where teeth are closely set.



SEPARATORS (UNCUT) 1,080/PKG.		
SEPARATOR TYPE	COLOR	.180"
SEP-O-LOOPS	ROYAL BLUE	500-001
	RADIO-OPAQUE (MIDNIGHT BLUE)	500-002
SEPO II	PURPLE (for widely spaced teeth)	500-003
	LIGHT BLUE (for closely spaced teeth)	500-004

RADIOOPAQUE SEPARATORS, LOOSE (CUT), BULK- 1,000/PKG.		
SIZE	COLOR	.180"
SMALL 0.040" Thick 1/8" O.D.	BLUE	500-005
	WHITE	500-006
LARGE 0.045" Thick 3/16" O.D.	BLUE	500-007
	WHITE	500-008

MINI-TIES

- .095 Diameter
- Stretchable Material
- Excellent Memory

- * Great for Single or Small Brackets - Holds Tight
- * Will Stretch around Various Size Brackets
- * Will Return Near Original Size - Holding Tight Around Bracket



Color	.095"
GRAY	600-001
CLEAR	600-002
SILVER	600-003

ELASTOMERIC CHAIN

Auradonics' chain elastics are designed with the same materials and properties as our elastomeric ligatures. This assures resistance to breakage and tearing, high flexibility and memory, color fastness, and the highest quality and consistency available to the orthodontic professional. Our three styles of chain spacing also provide more efficient and effective intraoral use. Chain colors match ligature tie colors. One 15 foot spool per package. Contain No Latex!

**Our chain colors match
the molded bracket ties !**



COLORS	SHORT	LONG	CONTINUOUS
GRAY	700-300	700-400	700-500
CLEAR	700-301	700-401	700-501
PINK	700-302	700-402	700-502
GREEN	700-303	700-403	700-503
PURPLE	700-304	700-404	700-504
RED	700-305	700-405	700-505
BLUE	700-306	700-406	700-506
BLACK	700-307	700-407	700-507
ORANGE	700-308	700-408	700-508
YELLOW	700-309	700-409	700-509
WHITE	700-310	700-410	700-510
IVORY	700-311	700-411	700-511
EMERALD	700-312	700-412	700-512
TURQUOISE	700-313	700-413	700-513
LIGHT BLUE	700-314	700-414	700-514
SILVER	700-315	700-415	700-515
GOLD	700-316	700-416	700-516
MAGENTA	700-317	700-417	700-517
KELLY GREEN	700-318	700-418	700-518
MARIGOLD	700-319	700-419	700-519
SPARKLE GREEN	700-320	700-420	700-520
SMOKE GRAY	700-321	700-421	700-521
TEAL	700-322	700-422	700-522
FUCHSIA	700-323	700-423	700-523
MIDNIGHT BLUE	700-324	700-424	700-524
GLOW IN DARK	700-325	700-425	700-525
PERIWINKLE	700-326	700-426	700-526
CAROLINA BLUE	700-327	700-427	700-527
BABY BLUE	700-328	700-428	700-528
SPARKLE RED	700-329	700-429	700-529
SPARKLE BLUE	700-330	700-430	700-530
BABY PINK	700-331	700-431	700-531
LILAC	700-332	700-432	700-532
PEARL	700-333	700-433	700-533
CORAL	700-334	700-434	700-534

FIGURE 8s - LIGATURES TO MAKE ORTHODONTISTS' WORK EASIER

One of our newest products is designed to make figure 8 ligating easier and better. Our "FIGURE 8" module innovation is designed to help even your novice assistants ligate an 8 without a tool in both hands. Works in horizontal and vertical orientation. Use vertical 8's to allow easier M-D movement (less friction!) and horizontal 8's when M-D is achieved to help

Figure 8s are sold in packs of 200, MEDIUM ONLY!

COLORS

GRAY	300-700	WHITE	300-710
CLEAR	300-701	TURQUOISE	300-713
NEON PINK	300-702	AQUA/LIGHT BLUE	300-714
NEON GREEN	300-703	SILVER	300-715
PURPLE	300-704	KELLY GREEN	300-718
RED	300-705	BABY BLUE	300-728
BLUE	300-706	BABY PINK	300-731
BLACK	300-707	LILAC	300-732
ORANGE	300-708	PEARL	300-733
YELLOW	300-709	CORAL	300-734

FIGURE 8 QUICK START
 Figure 8's come in 3 sizes, Small, Medium and Large. Choose the appropriate size for ligation based on the bracket size (Central, Lower Anterior, etc.) as well as how much force is required on that specific bracket. Make sure the Figure 8 is secured beneath tie-wings. Since brackets come in many shapes and sizes it is up to the practitioner to judge the fit is correct to remain in place until the next office visit.

Place ligature in this orientation to lessen frictional force on the archwire

Place ligature in this orientation to increase frictional force on the archwire



LIGATURES - LOADED ON CANES AND BULK

These latex-free ligature ties are pre-loaded on canes for doctors who find this dispensing system the most efficient. Also available in unloaded bulk quantities.

Our features include:

- ◆ A very high resistance to breakage or tearing during placement or treatment, saving valuable chair time and promoting patient confidence.
- ◆ High flexibility and resiliency without stiffness, making the ties easier to place and remove while providing a longer lasting, more effective tie. Our ties are colorfast to improve patient satisfaction and compliance.
- ◆ Available in the same variety of colors as our other ligatures.
- ◆ **Ask about bulk quantities of loose ligatures.**



COLORS

GRAY	300-400	KELLY GREEN	300-418
CLEAR	300-401	MARIGOLD	300-419
PINK	300-402	SPARKLE GREEN	300-420
NEON GREEN	300-403	SMOKE GRAY	300-421
PURPLE	300-404	TEAL	300-422
RED	300-405	FUCHSIA	300-423
ROYAL BLUE	300-406	MIDNIGHT BLUE	300-424
BLACK	300-407	GLOW IN DARK	300-425
ORANGE	300-408	PERIWINKLE	300-426
YELLOW	300-409	CAROLINA BLUE	300-427
WHITE	300-410	BABY BLUE	300-428
IVORY/TOOTH	300-411	SPARKLE RED	300-429
EMERALD	300-412	SPARKLE BLUE	300-430
TURQUOISE	300-413	BABY PINK	300-431
LIGHT BLUE	300-414	LILAC	300-432
SILVER	300-415	PEARL	300-433
GOLD	300-416	CORAL	300-434
MAGENTA	300-417		

